

# THG1L Antibody



PACO12739

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

THG1L

**Uniprot**

Q9NWX6

**Synonyms:**

tRNA-histidine guanylyltransferase 1-like (S. cerevisiae); THG1L; FLJ11601; FLJ20546; ICF45 ;

**Immunogen:**

Human THG1L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A