

# THNSL2 Antibody



PACO12740

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

THNSL2

**Uniprot**

Q86YJ6

**Synonyms:**

threonine synthase-like 2 (S. cerevisiae); THNSL2; FLJ10916; THS2; TSH2 ;

**Immunogen:**

Human THNSL2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A