

THOC7 Antibody



PACO12743

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

THOC7

Uniprot

Q6I9Y2

Synonyms:

THO complex 7 homolog (Drosophila); THOC7; FLJ23445; NIF3L1BP1; fSAP24 ;

Immunogen:

Human THOC7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A