

TIA1 Antibody



PACO12749

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TIA1

Uniprot

P31483

Synonyms:

TIA1 cytotoxic granule-associated RNA binding protein;

Immunogen:

Human TIA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A