

TIAF1 Antibody



PACO12750

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TIAF1

Uniprot

O95411

Synonyms:

TGFB1-induced anti-apoptotic factor 1; TIAF1; MAJN; MYO18A; MYSPDZ; SPR210 ;

Immunogen:

Human TIAF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A