

TICAM2 Antibody



PACO12751

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TICAM2

Uniprot

Q86XR7

Synonyms:

toll-like receptor adaptor molecule 2; TICAM2; MGC129876; MGC129877; TICAM-2; TIRAP3; TIRP; TRAM ;

Immunogen:

Human TICAM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A