

PACO12755

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TIMELESS

Uniprot

Q9UNS1

Synonyms:

timeless homolog (Drosophila); TIMELESS; FLJ12640; FLJ20714; TIM; TIM1; hTIM ;

Immunogen:

Human TIMELESS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A