

TIMM13 Antibody



PACO12757

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TIMM13

Uniprot

Q9Y5L4

Synonyms:

translocase of inner mitochondrial membrane 13 homolog (yeast); TIMM13; TIM13; TIM13B; TIMM13A; TIMM13B; ppv1 ;

Immunogen:

Human TIMM13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A