

TIMP2 Antibody



PACO12765

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TIMP2

Uniprot

P16035

Synonyms:

TIMP metalloproteinase inhibitor 2; TIMP2; CSC-21K ;

Immunogen:

Human TIMP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A