

# TLR4 Antibody



PACO12778

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat, Bovin

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TLR4

**Uniprot**

O00206

**Synonyms:**

toll-like receptor 4; TLR4; ARMD10; CD284; TOLL; hToll ;

**Immunogen:**

Human TLR4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A