

TM7SF3 Antibody



PACO12782

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TM7SF3

Uniprot

Q9NS93

Synonyms:

transmembrane 7 superfamily member 3;

Immunogen:

Human TM7SF3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A