

TMCC1 Antibody



PACO12784

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMCC1

Uniprot

O94876

Synonyms:

transmembrane and coiled-coil domain family 1; TMCC1; DKFZp686M0169; FLJ42680 ;

Immunogen:

Human TMCC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A