

# TMED3 Antibody



PACO12788

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMED3

**Uniprot**

Q9Y3Q3

**Synonyms:**

transmembrane emp24 protein transport domain containing 3; TMED3; C15orf22; MGC133022; P24B ;

**Immunogen:**

Human TMED3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A