

# TMEM11 Antibody



PACO12794

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMEM11

**Uniprot**

P17152

**Synonyms:**

transmembrane protein 11; TMEM11; C17orf35; PM1; PMI ;

**Immunogen:**

Human TMEM11.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A