

TMEM120A Antibody



PACO12795

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TMEM120A

Uniprot

Q9BXJ8

Synonyms:

transmembrane protein 120A; TMEM120A; TMPIT ;

Immunogen:

Human TMEM120A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A