

TMEM165 Antibody



PACO12797

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMEM165

Uniprot

Q9HC07

Synonyms:

transmembrane protein 165; TMEM165; TMPT27; TPRL ;

Immunogen:

Human TMEM165.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A