

# TMEM27 Antibody



PACO12802

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB, IHC, IF, ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMEM27

**Uniprot**

Q9HBJ8

**Synonyms:**

transmembrane protein 27; TMEM27; NX-17; NX17 ;

**Immunogen:**

Human TMEM27.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A