

# TMEM38A Antibody



PACO12803

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMEM38A

**Uniprot**

Q9H6F2

**Synonyms:**

transmembrane protein 38A; TMEM38A; MGC3169; TRIC-A; TRICA ;

**Immunogen:**

Human TMEM38A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A