

# TMEM50B Antibody



PACO12805

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMEM50B

**Uniprot**

P56557

**Synonyms:**

transmembrane protein 50B; TMEM50B; C21orf4; DKFZp686C2482; FLJ26146; HCVPT7TP3 ;

**Immunogen:**

Human TMEM50B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A