

TMEM67 Antibody



PACO12806

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

Protein Background:

Gene ID:

TMEM67

Uniprot

Q5HYA8

Synonyms:

112kd;

Immunogen:

Human TMEM67.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A