

TMLHE Antibody



PACO12808

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TMLHE

Uniprot

Q9NVH6

Synonyms:

trimethyllysine hydroxylase, epsilon; TMLHE; BBOX2; FLJ10727; TMLH; XAP130 ;

Immunogen:

Human TMLHE.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A