

TMPO Antibody



PACO12813

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TMPO

Uniprot

P42166

Synonyms:

thymopoietin; TMPO; CMD1T; LAP2; LEMD4; MGC61508; PRO0868; TP ;

Immunogen:

Human TMPO.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A