

TNFAIP6 Antibody



PACO12820

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TNFAIP6

Uniprot

P98066

Synonyms:

tumor necrosis factor, alpha-induced protein 6; TNFAIP6; TSG-6; TSG6 ;

Immunogen:

Human TNFAIP6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A