

# TNFAIP8L2 Antibody



PACO12822

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TNFAIP8L2

**Uniprot**

Q6P589

**Synonyms:**

tumor necrosis factor, alpha-induced protein 8-like 2; TNFAIP8L2; FLJ23467; TIPE2 ;

**Immunogen:**

Human TNFAIP8L2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A