

TNFRSF13B Antibody



PACO12824

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TNFRSF13B

Uniprot

O14836

Synonyms:

tumor necrosis factor receptor superfamily, member 13B; TNFRSF13B; CD267; CVID; FLJ39942; MGC133214; MGC39952; TACI; TNFRSF14B ;

Immunogen:

Human TNFRSF13B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A