

TNFSF12 Antibody



PACO12826

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TNFSF12

Uniprot

O43508

Synonyms:

tumor necrosis factor (ligand) superfamily, member 12; TNFSF12; APO3L; DR3LG; MGC129581; MGC20669; TWEAK ;

Immunogen:

Human TNFSF12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A