

# TNFSF13B Antibody



PACO12827

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TNFSF13B

**Uniprot**

Q9Y275

**Synonyms:**

tumor necrosis factor (ligand) superfamily, member 13b; TNFSF13B; BAFF; BLYS; CD257; DTL; TALL-1; TALL1; THANK; TNFSF20; ZTNF4 ;

**Immunogen:**

Human TNFSF13B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A