

TNFSF14 Antibody



PACO12828

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TNFSF14

Uniprot

O43557

Synonyms:

tumor necrosis factor (ligand) superfamily, member 14; TNFSF14; CD258; HVEM; LIGHT; LTg; TR2 ;

Immunogen:

Human TNFSF14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A