

TNIP2 Antibody



PACO12830

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TNIP2

Uniprot

Q8NFZ5

Synonyms:

TNFAIP3 interacting protein 2; TNIP2; ABIN-2; ABIN2; FLIP1; KLIP; MGC4289 ;

Immunogen:

Human TNIP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A