

TNPO2 Antibody



PACO12839

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TNPO2

Uniprot

O14787

Synonyms:

transportin 2; TNPO2; FLJ12155; IPO3; KPNB2B; TRN2 ;

Immunogen:

Human TNPO2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A