

TOR1A Antibody



PACO12854

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TOR1A

Uniprot

O14656

Synonyms:

torsin family 1, member A (torsin A); TOR1A; DQ2; DYT1 ;

Immunogen:

Human TOR1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A