

# TOR3A Antibody



PACO12856

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TOR3A

**Uniprot**

Q9H497

**Synonyms:**

torsin family 3, member A; TOR3A; ADIR; ADIR2; FLJ22345; MGC111104 ;

**Immunogen:**

Human TOR3A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A