

TP53 Antibody



PACO12857

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TP53

Uniprot

P04637

Synonyms:

tumor protein p53; TP53; FLJ92943; LFS1; TRP53; p53 ;

Immunogen:

Human TP53.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A