

TPD52L1 Antibody



PACO12860

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TPD52L1

Uniprot

Q16890

Synonyms:

tumor protein D52-like 1; TPD52L1; D53; MGC8556; TPD52L2; hD53 ;

Immunogen:

Human TPD52L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A