

TPT1 Antibody



PACO12877

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TPT1

Uniprot

P13693

Synonyms:

tumor protein, translationally-controlled 1; TPT1; FLJ27337; HRF; TCTP; p02 ; Fortilin; p23

Immunogen:

Human TPT1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A