

TRA2A Antibody



PACO12880

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TRA2A

Uniprot

Q13595

Synonyms:

transformer 2 alpha homolog (Drosophila); TRA2A; HSU53209 ;

Immunogen:

Human TRA2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A