

TRAF1 Antibody



PACO12882

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRAF1

Uniprot

Q13077

Synonyms:

TNF receptor-associated factor 1; TRAF1; EBI6; MGC:10353 ;

Immunogen:

Human TRAF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A