

TRAF5 Antibody



PACO12885

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRAF5

Uniprot

O00463

Synonyms:

TNF receptor-associated factor 5; TRAF5; MGC:39780; RNF84 ;

Immunogen:

Human TRAF5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A