

# TRAF6 Antibody



PACO12886

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRAF6

**Uniprot**

Q9Y4K3

**Synonyms:**

TNF receptor-associated factor 6; TRAF6; MGC:3310; RNF85 ;

**Immunogen:**

Human TRAF6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A