

# TRAK2 Antibody



PACO12889

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRAK2

**Uniprot**

O60296

**Synonyms:**

trafficking protein, kinesin binding 2; TRAK2; ALS2CR3; CALS-C; GRIF-1; GRIF1; KIAA0549; OIP98 ;

**Immunogen:**

Human TRAK2.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A