

TRAM2 Antibody



PACO12891

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRAM2

Uniprot

Q15035

Synonyms:

translocation associated membrane protein 2; TRAM2; KIAA0057 ;

Immunogen:

Human TRAM2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A