

TRAPPC9 Antibody



PACO12896

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRAPPC9

Uniprot

Q96Q05

Synonyms:

trafficking protein particle complex 9; TRAPPC9; IBP; IKBKBBP; KIAA1882; MGC4737; MGC4769; NIBP; T1; TRS120 ;

Immunogen:

Human TRAPPC9.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A