

TRAT1 Antibody



PACO12897

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRAT1

Uniprot

Q6PIZ9

Synonyms:

T cell receptor associated transmembrane adaptor 1; TRAT1; HSPC062; TCRIM; TRIM ;

Immunogen:

Human TRAT1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A