

TRIM14 Antibody



PACO12903

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TRIM14

Uniprot

Q14142

Synonyms:

tripartite motif-containing 14; TRIM14; KIAA0129 ;

Immunogen:

Human TRIM14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A