

# TRIM23 Antibody



PACO12908

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

TRIM23

**Uniprot**

P36406

**Synonyms:**

tripartite motif-containing 23; TRIM23; ARD1; ARFD1; RNF46 ;

**Immunogen:**

Human TRIM23.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A