

TRIM37 Antibody



PACO12915

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRIM37

Uniprot

O94972

Synonyms:

tripartite motif-containing 37; TRIM37; KIAA0898; MUL; POB1; TEF3 ;

Immunogen:

Human TRIM37.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A