

TRIM39 Antibody



PACO12917

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRIM39

Uniprot

Q9HCM9

Synonyms:

tripartite motif-containing 39; TRIM39; MGC32984; RNF23; TFP; TRIM39B ;

Immunogen:

Human TRIM39.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A