

TRIM54 Antibody



PACO12922

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRIM54

Uniprot

Q9BYV2

Synonyms:

tripartite motif-containing 54; TRIM54; MURF; MURF-3; RNF30 ;

Immunogen:

Human TRIM54.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A