

TRIM69 Antibody



PACO12924

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TRIM69

Uniprot

Q86WT6

Synonyms:

tripartite motif-containing 69; TRIM69; HSD34; RNF36; Trif ;

Immunogen:

Human TRIM69.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A