

TRIOBP Antibody



PACO12926

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TRIOBP

Uniprot

Q9H2D6

Synonyms:

TRIO and F-actin binding protein; TRIOBP; DFNB28; FLJ39315; HRIHFB2122; KIAA1662; TARA; dJ37E16.4 ;

Immunogen:

Human TRIOBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A